

# monthly 07

LIMITED TIME OFFERS! July

## Midway Select

FREE product automatically ships with order!

DOUBLE  
CONNECT  
CASH



**MIDWAY**  
SELECT

**Midway Select Temporary Cement:** Non-Eugenol, residue free, superb adhesion easy removal, neutral color

**002-0648** Midway Select Temporary Cement NE 5mL automix syringe + 10 mixing tips

**Buy 2 Temporary Cement**

**002-0580** Midway Select HP Mixing Tips Brown Pkg/20

**Get 1 Mixing tip Free**

**Midway Select Flowable Composite:** Universal nanohybrid, low viscosity, light cure restorative

**002-0584** Midway Everyday Flowable Composite A2 4 x 1.5g syringes + 20 tips

**Buy 3 flowable, mix & match**

**002-0585** Midway Everyday Flowable Composite A3 4 x 1.5g syringes + 20 tips

**002-0586** Midway Everyday Flowable Composite A3.5 4 x 1.5g syringes + 20 tips

**002-0587** Midway Everyday Flowable Composite B1 4 x 1.5g syringes + 20 tips

**002-0608** Midway Select Etch 40% – 30 ml Bulk Syringe with Docking Port

**Get 1 Etch Free**

DOUBLE  
CONNECT  
CASH



## NEW Product Feature



### PANAVIA Veneer LC Cement

PANAVIA™ Veneer LC is an advanced light-cure resin cement with excellent color stability, high wear resistance, provides ample working time, and offers fast and easy cleanup. It is indicated for permanent cementation of aesthetic veneers and can also be used for inlays and onlays less than 2mm thick. Offering an ideal consistency for easy and secure application across the intaglio surface, PANAVIA™ Veneer LC does not stick to the application tip due to its unique nanocluster filler technology. This technology consists of spherical silica fillers in the cement to provide stability and control, while simultaneously promoting flow when the veneer is seated on the tooth, thereby allowing easy veneer placement without drifting or sagging. PANAVIA™ Veneer LC is available in four highly aesthetic shades, all featuring great polishability and an excellent gloss retention to support aesthetic margins over time. Matching try-in pastes (PANAVIA™ V5 Try-In Paste) allow for a simulation of the final appearance prior to final cementation.

**906-15167** Panavia Veneer LC Cement Intro Kit

**906-15169** Panavia Veneer LC Cement 2.6g Syringe - Clear

**906-15171** Panavia Veneer LC Cement 2.6g Syringe - White

**906-15168** Panavia Veneer LC Cement 2.6g Syringe - Universal A2

**906-15170** Panavia Veneer LC Cement 2.6g Syringe - Brown A4

## Connect Cash Bonus



### Air Techniques Disposable Nasal Hoods (Box of 12)

Buy 1 box,  
Get 1 box of the same FREE!

#### Adult

- 906-14401 Blueberry Dance
- 906-14402 Strawberry Delight
- 906-14406 Non-Scented Pure Ghost
- 906-14405 Chilla Vanilla

#### Pediatric

- 906-14409 Bouncy Bubble
- 906-14407 Blueberry Dance
- 906-14408 Strawberry Delight
- 906-14412 Non-Scented Pure Ghost
- 906-14411 Chilla Vanilla

Promo Code **NitrousSummer22**

## Equipment of the Month

### Enbio Sterilizer



Why would you waste an hour on something you can do in 7 minutes? Enbio is the fastest sterilizer on the market. It is also gloriously intuitive and easy to use. It will fit any interior thanks to its compact size (25 x 16 x 56 cm). At the same time, it has the space for all the standard tools, in their packaging or without it. Fancy a chat standing next to the running Enbio? No problem. It is so quiet you can relax by it (32dB, comparable to an empty, tranquil road).

**\$2,899.00**

DOUBLE  
CONNECT  
CASH



## Midway EXCLUSIVE Promotions



### Sheer White

- 906-14875 Sheer White Take Home Whitening Kit, 6 Patient/pkg
- 906-14876 Sheer White Take Home Whitening Kit, 24 Patient/pkg

**Buy 1 Get 1 of  
the same FREE!**



#### DentaPure Municipally Plumbed Dental Waterline Treatment Cartridge - Crosstex

Annual installation in dental units directly plumbed to city water (via included poly-flo fittings) provides compliant dental procedural water meeting current EPA standard for potable water.

901-5175

**Buy 2 Get 1 FREE!**

FREE product automatically ships with order!

#### Specifications:

- 1 installation provides 365 days (240L) of procedural water than meets the current EPA standard for potable water
- Minimal maintenance following installation
- Cartridge reduces patient and staff bacterial exposure by eluting non-allergenic elemental iodine into dental treatment water
- Change every 365 days or 240L of water
- Use NOT restricted by the EPA Rule BMP for amalgam waste

Manufacturer ID: DP365M10



#### DentaPure Bottle Waterline Treatment Cartridge - Crosstex

Annual installation onto the dental unit water bottle pick-up tube via included luer fitting provides compliant dental procedural water meeting current EPA standard for potable water.

389-1500

**Buy 2 Get 1 FREE!**

FREE product automatically ships with order!

#### Specifications:

- 1 installation provides 365 days (240L) of procedural water than meets the current EPA standard for potable water
- For use for potable water
- Quickly and easily installed by office staff
- Minimal maintenance following installation
- Cartridge reduces patient and staff bacterial exposure by eluting non-allergenic elemental iodine into dental treatment water
- Change every 365 days or 240L of water
- Use NOT restricted by the EPA Rule BMP for amalgam waste

Manufacturer ID: DP365B10

**HuFriedyGroup**  
The Best In Practice